

QuickGuide: RealFast™ Variant Detection on Roche LightCycler® 480II

Setup for Variant Detection and Onco RealFast™ Assays:

- Open the **LightCycler® 480 software** and login with your username and password.
 - Choose **New Experiment**.
- Define your PCR program in the **Run Protocol** tab.
 - Select **Dual Color Hydrolysis Probe / UPL Probe** as Detection Format.
 - Select a **reaction volume** of 20µl and setup the **Program**:
 - Add programs in the **Program Name** window with „+“ and edit **Cycles** and **Analysis Mode**. To edit **Acquisition Mode** and **Hold**, click on the corresponding step in the **Program Name** window and change parameters or add steps („+“) in the **Temperature Targets** window.
 - Setup the PCR program according to the Instructions for Use of your RealFast™ Assay.
- Click on **Subset Editor** button on the left side of the window.
 - Press „+“ to create a new subset and rename your subset.
 - Select wells in the grid and press **Apply**.





- Click on **Sample Editor** button on the left side of the window.

- **Select Workflow > Abs Quant.**
- **Select Filter Combinations > 456-510nm (FAM) and 533-580nm (HEX).**
- Choose your **Subset** of Samples.







- Define your **No Template Control (NTC)**:

- **Select Samples** field: select well by mouse-click or two wells by ctrl+mouse click.
- **Edit Abs Quant Properties** field: Type **NTC** in the **Sample Name** field and press **Enter**. Choose **Negative Control** as **Sample Type**.

Pos	Filter combination	Color	Repl Of	Sample Name	Quantification Sample Type	Concentration
▶ A1	FAM (465-510)		A1	NTC	Negative C	
A1	VIC / HEX /		A1	NTC	Negative Con	
A2	FAM (465-510)		A1	NTC	Negative Con	
A2	VIC / HEX /		A1	NTC	Negative Con	

- Define your **Samples**:

- **Select Samples** field: select a well by mouse click.
- **Edit Abs Quant Properties** field: Type the name of your first sample in the corresponding field and press **Enter**. Check **Unknown**.
- Define the rest of your samples.

Pos	Filter combination	Color	Repl Of	Sample Name	Quantification Sample Type	Concentration
▶ A3	FAM (465-510)		A3	3	Unknown	
A3	VIC / HEX /		A3	3	Unknown	
A4	FAM (465-510)		A3	3	Unknown	
A4	VIC / HEX /		A3	3	Unknown	


- **Save** or **export** your experiment by pressing the corresponding button:





- Load your samples and start the experiment.




Analysis of Variant Detection and Onco RealFast™ Assays using “Abs Quant/Fit Points”:

- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
 - Choose an experiment from the data bank or
 - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.

- Within the **Create New Analysis** field choose **Abs Quant/Fit Points**.
- A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button
- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Within the **Noise Band** field adjust the threshold manually. Set the threshold in the exponential phase of the amplification curves above the background.





Abs Quant results				
Positive Negative Standard				
Samples				
Include	Color	Pos	Name	Cp
<input checked="" type="checkbox"/>	Green	A4	B*070201	
<input checked="" type="checkbox"/>	Green	B4	B*070201	
<input checked="" type="checkbox"/>	Red	C4	B*2701	24, 45
<input checked="" type="checkbox"/>	Red	D4	B*2701	24, 59
<input checked="" type="checkbox"/>	Green	E4	H2O	
<input checked="" type="checkbox"/>	Green	F4	H2O	




- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.
- To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
- Press the **Calculate** button and review your results in the chart bottom-left. All genomic DNA samples and Controls show an amplification in the HEX channel and are therefore marked as positive.


» **Note:** For **Onco RealFast™** Assays (e.g. **EGFR T790M**,...), it is **mandatory** to do the analysis according to the section **Data Analysis / Interpretation of Results** in the respective **Instructions for Use**. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.

Optional:

Analysis of Variant Detection and Onco RealFast™ Assays using Abs Quant/2nd Derivative Max:


- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
 - Choose an experiment from the data bank or
 - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.
 - Within the **Create New Analysis** field choose **Abs Quant/2nd Derivative Max**.
 - A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button.




Create new analysis

Analysis Type	* Abs Quant/2nd Derivative Max
Subset	* HLA B57
Program	* amplification
Name	* Abs Quant/2nd Derivative Max for HLA B57



- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.



Abs Quant results				
Legend: Positive (red), Negative (green), Uncertain (blue), Standard (dark red)				
Samples				
Include	Co...	Pos	Name	Cp
<input checked="" type="checkbox"/>	█		A1 NTC	
<input checked="" type="checkbox"/>	█		A2 NTC	
<input checked="" type="checkbox"/>	█		A5 4	
<input checked="" type="checkbox"/>	█		A6 4	
<input checked="" type="checkbox"/>	█	█	A9 9	25, 89
<input checked="" type="checkbox"/>	█	█	A10 9	25, 81



- To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
 - Press the **Calculate** button and review your results in the chart bottom-left. All DNA samples and the Positive Control show an amplification in the HEX channel and are therefore marked as positive.
- » **Note:** For **Onco RealFast™ Assays** (e.g. EGFR T790M,...), it is **mandatory** to do the analysis according to the section **Data Analysis / Interpretation of Results** in the respective **Instructions for Use**. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.