


## QuickGuide: RealFast™ Variant Detection on LightCycler® 480II

### Setup for Variant Detection and Onco RealFast™ Assays:

- Open the **LightCycler® 480 software** (QuickGuide is based on version 1.5.1) and login with your username and password.
  - Choose **New Experiment**.
- Define your PCR program in the **Run Protocol** tab.
  - Select **Dual Color Hydrolysis Probe / UPL Probe** as Detection Format.
  - Select a **reaction volume** of 20 µl and setup the **Program**:
  - Add programs in the **Program Name** window with "+" and edit **Cycles** and **Analysis Mode**. To edit **Acquisition Mode** and **Hold**, click on the corresponding step in the **Program Name** window and change parameters or add steps ("+") in the **Temperature Targets** window.
  - Setup the PCR program according to the Instructions for Use of your RealFast™ Assay.
- Click on **Subset Editor** button on the left side of the window.
  - Press "+" to create a new subset and rename your subset.
  - Select wells in the grid and press **Apply**.
- Click on **Sample Editor** button on the left side of the window.
  - **Select Workflow > Abs Quant**.
  - **Select Filter Combinations > 456-510 nm (FAM) and 533-580 nm (HEX)**.
  - Choose your **Subset** of Samples.



- Define your **No Template Control (NTC)**:

- **Select Samples** field: select well by mouse-click or two wells by ctrl + mouse click.
- **Edit Abs Quant Properties** field: Type **NTC** in the **Sample Name** field and press **Enter**. Choose **Negative Control** as **Sample Type**.

Pos	Filter combination	Color	Repl Of	Sample Name	Quantification Sample Type	Concentration
▶ A1	FAM (465-510)	<span style="color: blue;">■</span>	A1	NTC	Negative C	
A1	VIC / HEX /	<span style="color: blue;">■</span>	A1	NTC	Negative Con	
A2	FAM (465-510)	<span style="color: red;">■</span>	A1	NTC	Negative Con	
A2	VIC / HEX /	<span style="color: red;">■</span>	A1	NTC	Negative Con	

- Define your **Samples**:

- **Select Samples** field: select a well by mouse click.
- **Edit Abs Quant Properties** field: Type the name of your first sample in the corresponding field and press **Enter**. Check **Unknown**.
- Define the rest of your samples.

Pos	Filter combination	Color	Repl Of	Sample Name	Quantification Sample Type	Concentration
▶ A3	FAM (465-510)	<span style="color: green;">■</span>	A3	3	Unknown	
A3	VIC / HEX /	<span style="color: green;">■</span>	A3	3	Unknown	
A4	FAM (465-510)	<span style="color: magenta;">■</span>	A3	3	Unknown	
A4	VIC / HEX /	<span style="color: magenta;">■</span>	A3	3	Unknown	


- **Save** or **export** your experiment by pressing the corresponding button:



- Load your samples and start the experiment.




## Analysis of Variant Detection and Onco RealFast™ Assays using “Abs Quant/Fit Points”:

- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
  - Choose an experiment from the data bank or
  - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.



- Within the **Create New Analysis** field choose **Abs Quant/Fit Points**.
- A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button

- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Within the **Noise Band** field adjust the threshold manually. Set the threshold in the exponential phase of the amplification curves above the background.

Abs Quant results				
<input type="checkbox"/> Positive <input type="checkbox"/> Negative <input type="checkbox"/> Standard				
Samples				
Include	Color	Pos	Name	Cp
<input checked="" type="checkbox"/>	<span style="color: green;">■</span>	A4	B*070201	
<input checked="" type="checkbox"/>	<span style="color: green;">■</span>	B4	B*070201	
<input checked="" type="checkbox"/>	<span style="color: red;">■</span>	C4	B*2701	24, 45
<input checked="" type="checkbox"/>	<span style="color: red;">■</span>	D4	B*2701	24, 59
<input checked="" type="checkbox"/>	<span style="color: green;">■</span>	E4	H2O	
<input checked="" type="checkbox"/>	<span style="color: green;">■</span>	F4	H2O	





- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.
  - To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
  - Press the **Calculate** button and review your results in the chart bottom-left. All genomic DNA samples and Controls show an amplification in the HEX channel and are therefore marked as positive.
- » **Note:** Set threshold value for the FAM channel just above the background fluorescent signal generated by the e.g. HLA-B27 Negative Control. For **Onco RealFast™ Assays** (e.g. EGFR T790M), it is **mandatory** to do the analysis according to the section *Data Analysis / Interpretation of Results* in the respective *Instructions for Use*. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.


Optional:

### **Analysis of Variant Detection and Onco RealFast™ Assays using Abs Quant/2nd Derivative Max:**


- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
  - Choose an experiment from the data bank or
  - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.
  - Within the **Create New Analysis** field choose **Abs Quant/2nd Derivative Max**.
  - A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button.





- *Optional:* Press **Color Compensation** and choose **In Use** for color compensation of FAM (510) and VIC (580). Within the pop-up window select a color compensation and press the **OK** button.
- Press **Filter Comb 465-510**, choose **FAM (465-510)** in the pop-up window and press the **OK** button.
- Press the **Calculate** button. Choose **Abs Quant Results** in the chart bottom-left and review your samples. Samples positive for the targeted variant show an amplification in the FAM channel and will be color-coded.



Include	Co...	Pos	Name	Cp
<input checked="" type="checkbox"/>		Green	A1 NTC	
<input checked="" type="checkbox"/>		Green	A2 NTC	
<input checked="" type="checkbox"/>		Green	A5 4	
<input checked="" type="checkbox"/>		Green	A6 4	
<input checked="" type="checkbox"/>		Red	A9 9	25,89
<input checked="" type="checkbox"/>		Red	A10 9	25,81



- To review proper amplification of all samples, press **Filter Comb 465-510**, choose **VIC / HEX / Yellow 555 (533/580)** in the pop-up window and press the **OK** button.
  - Press the **Calculate** button and review your results in the chart bottom-left. All DNA samples and the Positive Control show an amplification in the HEX channel and are therefore marked as positive.
- » **Note:** Set threshold value for the FAM channel just above the background fluorescent signal generated by the e.g. HLA-B27 Negative Control. For **Onco RealFast™ Assays** (e.g. EGFR T790M), it is **mandatory** to do the analysis according to the section *Data Analysis / Interpretation of Results* in the respective *Instructions for Use*. «

- **Save** or **export** your data by pressing the corresponding button.  or 
- After saving your data you can customize and generate a report via the **Report** button.