

QuickGuide: RealFast™ Genotyping on AB QuantStudio 5

Setup for Genotyping Assays:


- Open the QuantStudio Design & Analysis Software (QuickGuide is based on version 1.5.2) and click **Create New Experiment**.

- Define **Experiment Properties**
 - Name: Name of experiment
 - Instrument type: **QuantStudio 5 System**
 - Block type: according to your instrument
 - Experiment type: **Genotyping**
 - Chemistry: **TaqMan® Reagents**
 - Run mode: **Standard**

- In **Method** select a sample volume of **20 µl** and setup the PCR program:
optional but recommended: include Pre-PCR Read 30 sec at 60°C.

Include: Pre-PCR Read Amplification Post-PCR Read



- Holding Stage: **3 min at 95°C**
- Cycling Stage: **40 cycles of 15 sec at 95°C and 1 min at 60°C.**
Make sure **Data Collection On** is enabled.



- In **Plate Quick Setup** define **None** as **Passive Reference** dye.




- In **Plate Advanced Setup** assign **SNP Assay(s)** and **Samples** to selected wells.
 - Select **SNP Assay 1** and click **Action** and **Edit** and enter assay name. Keep default settings for reporter dyes:
 - Allele 1 = VIC (corresponds to HEX-labeled probe)
 - Allele 2 = FAM (corresponds to FAM-labeled probe)
 - Quencher = NFQ-MGB




» **Note:** In case you want to run several **different singleplex RealFast™ Assays in parallel**, create **New SNP Assay** for each of the markers to be analyzed. «

» **Note:** For **multiplex RealFast™ Assays two SNP Assays have to be created** – one for each marker. «



- Click **Add** to open a **New SNP Assay** and enter assay name. Change reporter dyes:
 - Allele 1 = CY5 (corresponds to CY5-labeled probe)
 - Allele 2 = ROX (corresponds to ROX-labeled probe)
 - Quencher = NFQ-MGB



- Click **Add** repeatedly to enter all your samples and controls.



- Select the total number of wells per assay (2 Negative Controls + Positive Controls + number of samples) by click+drag in the grid in the Plate Layout.
- Assign SNP assay to selected wells by ticking the box next to the **SNP Assay(s)**.

» **Note:** For **multiplex RealFast™** Assays **assign both SNP Assays** to the selected wells. «



- Define your **Negative Control**:

- Select a replicate (2 wells) in the plate layout by click+drag.
- Assign Negative Control to corresponding well by ticking the box next to the sample name.
- Select "N" in the pull-down menu **Task**.



- *Optional:* Define your **Positive Controls**:

- In the plate layout select a well for each **Positive Control** by mouse click.
- Assign Positive Control to corresponding well by ticking the box next to the sample name.
- In the pull-down menu **Task** select
 - “1/1” for **HEX-** or **Cy5-positive Control** (in most cases WT-Control),
 - “2/2” for **FAM-** or **ROX-positive Control** (in most cases MUT-Control),
 - “1/2” for **HEX-/FAM-** or **CY5-/ROX-positive Control** (mix WT- and MUT-Control 1:1)

Alternatively, Positive controls can be defined as unknown (“U”) (see below).

» **Note:** For assignment of fluorophores to hydrolysis probes see *Instructions for Use of the respective RealFast™ Assay*. «

- Define your **Samples**:

- In the plate layout select a well for each sample by mouse click.
- Assign sample to corresponding well by ticking the box next to the sample name.
- Select “U” in the **Task** pull-down menu.

- Load your reaction plate into the QuantStudio 5 instrument. Go to **Run** and press **START RUN**, select the serial number of the machine to start the run.

Analysis of Genotyping Assays:

After completing a run or after opening an **Existing Experiment** the software displays the **Results**:

- The **Amplification Plot** appears automatically.
- If you have **different RealFast™** Assays in one run, select appropriate samples by clicking **Action** and selecting **Select SNP Assay** in the dropdown menu.


» **Note:** With **multiplex RealFast™** Assays all 4 fluorophores are depicted in one **Amplification Plot!**

- To change the threshold according to the IFU click **All Targets**, select a target and move the threshold with the cursor, then click **Analyze**.




- **Allelic Discrimination Plot:** Cartesian Plot type


- Always verify correct assignment of samples in the **Allelic Discrimination Plot:** e.g. select WT-Control in the **Plate Advanced Setup** > corresponding point on the Allele 1 axis turns from red to turquoise.
- In case a sample appears as undetermined = x, verify correct amplification in the **Amplification Plot** and manually assign genotype with **Apply Call** in dropdown menu at the eye icon.



- Click **Amplification Plot** to control correct amplification of all controls:

- Select the following plot settings at the eye icon:
 - > **Plot Type:** ΔRn vs Cycle
 - > **Graph Type:** Linear
 - > **Plot Color:** Target.
- In the dropdown menu **All Target** choose **Allele 1**. In **View Plate Layout** select your **Positive Control** for **Allele 1** (mostly WT-Control) - an amplification curve should be visible in the Amplification Plot. Select your **Positive Control** for **Allele 2** (mostly MUT-Control) - NO amplification curve should be visible.
- The threshold for **Allele 1** should be above the background signal of the **Positive Control** for **Allele 2**. If not, click on the eye icon and disable tickbox for **Auto** Threshold. Set the threshold manually by clicking on the threshold line in the plot and moving it above the background signal.
- In the dropdown menu **All Target** choose **Allele 2**. Select your **positive Control** for **Allele 2** (mostly MUT- Control) - an amplification curve should be visible in the Amplification Plot. Select your **Positive Control** for **Allele 1** (mostly WT-Control) - NO amplification curve should be visible.
- The threshold for **Allele 2** should be above the background signal of the **Positive Control** for **Allele 1**. If not proceed as described above for threshold setting.
- Verify absence of any contamination in the **Negative Control**. No amplification should be visible, neither for Allele 1 nor for Allele 2.
- Choose **All Target**. Select your samples one by one and verify positive amplification.



- To show results as table: click 

- Adjust the table according to your needs by selecting/deselecting the listed features in **View**.

- To print a report select **File > Print Report** in the upper menu bar:

- Select data for the report according to your needs.
- If you wish to export your data go to the **Export** tab and choose the file type and content