



# StripAssay<sup>®</sup> Technology & Procedure

## **ViennaLab StripAssays**<sup>®</sup>

Human diagnostic assays for genotyping mutations and polymorphisms (SNPs, deletions, insertions) based on **polymerase chain reaction (PCR)** and reverse-hybridization to allele-specific oligonucleotide probes immobilized on **teststrips**.

## StripAssay® components



SIMPLE

ECONOMICAL

RELIABLE

RAPID

EFFICIENT

CONVENIENT

Ready-to-use Amplification Mixes, Teststrips, processing reagents and incubation trays.

## StripAssay<sup>®</sup> components

### **Materials required but not supplied\*):**

- adjustable microcentrifuge capable of 3,000-12,000 rpm (1,000-12,000 x g)
- incubator (e.g. heating block, water bath) capable of 56°C and 98°C ( $\pm 2^\circ\text{C}$ )
- thermocycler and suitable thin-walled plastic reaction tubes/strips
- waterbath with shaking platform and adjustable temperature ( $45^\circ\text{C} \pm 1^\circ\text{C}$ )  
alternatively: thermoshaker with heated lid (selected models only!) plus simple non-shaking waterbath
- vacuum aspiration apparatus
- shaker (rocker or orbital shaker)
- for most StripAssays<sup>®</sup>: Taq DNA polymerase (avoid hot-start enzymes!)
- for DNA from external kits/protocols: spectrophotometer (UV 260 nm)  
for DNA extracted from FFPE tissue: fluorometer (e.g. Qubit system)
- optional: agarose gel electrophoresis equipment  
(for control of amplification products)

*\*) In addition to standard molecular biology laboratory equipment*

## StripAssay<sup>®</sup> procedure

1. **DNA Isolation**
2. ***In Vitro* Amplification (PCR)**
3. **Hybridization**
4. **Stringent Washes**
5. **Color Development**
6. **Interpretation**


# 1. DNA Isolation from blood (GEN<sup>X</sup>TRACT)

- Use **fresh** or **frozen blood** (anticoagulant: EDTA; avoid citrate or heparin)  
max. **3 days** at **ambient temperature**  
max. **1 week** at **2-8°C**  
max. **1 year** **frozen**
- **Avoid** repeated **freeze-thaw** cycles
- Bring blood samples to **room temperature**
- Before withdrawing an aliquot, **mix well** by inverting blood collection tubes several times
- Allow Lysis Solution and GEN<sup>X</sup>TRACT Resin to reach **room temperature**

Pipette **100  $\mu$ l blood sample** into a 1.5 ml microtube with screw cap.



Add **1 ml Lysis Solution**, close tube and mix by inverting several times.



Let stand for **15 min**  
at room temperature.





Centrifuge for 5 min at 3,000 rpm (approx. 1,000 x g) in a micro-centrifuge.




Remove and discard the upper (top) 1 ml of supernatant.

*Leave pellet undisturbed !*





Add **1 ml Lysis Solution**, close tube and mix by inverting several times.



*Pellet should be completely resuspended !*

Centrifuge for **5 min**  
at **12,000 rpm**  
(approx. 12,000 x g)  
in a microcentrifuge.



Remove and discard the supernatant except for approx. **50 µl of a visible, soft pellet.**



Resuspend GEN<sup>X</sup>TRACT Resin  
by swirling the bottle thoroughly.

*GEN<sup>X</sup>TRACT Resin sediments quickly.  
Repeat resuspension each time immediately  
before removing another aliquot !*



Add **200  $\mu$ l GEN<sup>X</sup>TRACT Resin** to the pellet.



Close tube and vortex for 10 sec.



**Incubate for 20  
min at 56°C.**

Vortex for 10 sec.



Incubate for **10 min**  
at **98°C**.

Vortex for 10 sec.





Centrifuge for **5 min** at **12,000 rpm** in a microcentrifuge. Cool on ice.



*The resulting supernatant contains DNA template suitable for immediate use in PCR.*

*For further storage, the supernatant should be transferred into a fresh tube and kept refrigerated (2-8°C; up to one week) or frozen.*



## 2. *In Vitro* Amplification (PCR)

- Use **10-50 ng** DNA per reaction (**2-10 ng/μl**), if DNA templates are not prepared by the GEN<sup>x</sup>TRACT protocol  
(DNA from FFPE tissue: use **1-10 ng/μl** quantified by fluorometry)
- Keep all **PCR reagents** and **DNA templates refrigerated**
- Prepare one reaction tube for each sample and place tubes on ice
- Perform all steps **on ice** (0-4°C) until start of the thermocycler
- **Preheat** thermocycler to starting temperature (usually 94°C)

Prepare a fresh working dilution (usually 0.2 U/ $\mu$ l) of **Taq DNA polymerase** in **Taq Dilution Buffer** (transparent cap).



For each sample prepare a final PCR reaction mix on ice:

**15  $\mu$ l Amplification Mix** (yellow or green or white cap)



**+ 5  $\mu$ l diluted Taq DNA polymerase** (usually 1U)



**+ 5  $\mu$ l DNA template**



Cap tubes tightly.  
Insert reaction tubes into the  
preheated thermocycler and  
run the appropriate program.

*Store amplification products on ice  
or at 2-8°C for further use.*



## 3. Hybridization (at 45°C)

- Use a **shaking waterbath**
- Adjust water level to approx. **½ of the height of the Typing Tray**
- Set temperature to exactly **45°C**
- Check water temperature with a fully inserted calibrated thermometer



..... or alternatively:

Use **thermoshaker with heated lid** (for incubations) plus a simple waterbath (for preheating reagents).




Set temperature to exactly 45°C.

**Do not use air-based hybridization devices because of insufficient heat transmission (e.g. hybridization oven).**





*Prewarm Hybridization Buffer and Wash Solution A to 45°C.*




*Take care that all precipitates formed at 2-8°C become completely dissolved after prewarming.*



*Allow Teststrips, DNAT, Conjugate Solution, Wash Solution B and Color Developer to reach room temperature. Prepare Typing Tray(s).*

*Remove one Teststrip for each sample using clean tweezers. (Touch Teststrips with gloves only!)*



*Label Teststrips outside of the marker lines with a pencil. (No ballpoint pens, ink markers, etc.)*

Pipette **10 µl DNAT** (blue cap) into the lower corner of each lane to be used in the Typing Trays (one lane per sample).

Add **10 µl amplification product** into the corresponding drop of DNAT. Mix thoroughly with a pipette.

*(The solution will remain blue.)*

Let stand for **5 min** at room temperature.



**Add 1 ml Hybridization Buffer**  
(prewarmed to 45°C) into each  
lane. Gently agitate tray.

*The blue color will disappear !*



Insert **Teststrips** with marked side up (lines visible!) into the respective lanes.


Submerge Teststrips completely.

*In order to avoid cross-contamination, do **not** dip tweezers into solution.*



Incubate for **30 min** at **45°C** on the shaking platform of the waterbath or thermoshaker.

*Set moderate shaking frequency to avoid spilling. Keep the cover of the waterbath or thermoshaker closed to avoid variations in temperature.*



At the end of incubation  
remove hybridization  
solutions by vacuum  
aspiration.

*Proceed immediately.  
Do not allow Teststrips to run  
dry during the entire procedure.*



## 4. Stringent Washes (at 45°C)

One **rinsing** step:

Add **1 ml Wash Solution A** prewarmed to 45°C. Rinse briefly (10 sec). Remove liquids by vacuum aspiration.

Two **washing** steps:

Add **1 ml Wash Solution A** (45°C).  
Incubate for **15 min** at **45°C** on the shaking platform of the waterbath or thermoshaker.  
Remove liquids by vacuum aspiration.

**Repeat 1 x**



## 5. Color Development (at room temperature)

Add **1 ml Conjugate Solution**.  
Incubate for **15 min at room temperature** on a rocker or orbital shaker. Remove liquids by vacuum aspiration.



One **rinsing** step:

Add **1 ml Wash Solution B**. Rinse briefly (10 sec). Remove liquids by vacuum aspiration.


Two **washing** steps:

Add **1 ml Wash Solution B**. Incubate for 5 min at room temperature on a rocker or orbital shaker. Remove liquids by vacuum aspiration.

**Repeat 1 x**



**Add 1 ml Color Developer.**  
**Incubate for 15 min at room temperature**  
**in the dark** on a rocker or orbital shaker.




*A purple staining will appear upon positive reaction.*



Wash Teststrips several times with distilled water.



Let strips dry in the dark on absorbent paper.



*Do not expose Teststrips to intense light after Color Development.*

- 3. Hybridization
- 4. Stringent Washes
- 5. Color Development

## *Optional:* Automated processing of teststrips

Three steps on one instrument: Hybridization + Stringent Washes + Color Development fully automated under controlled temperature

for 48 teststrips:

**ProfiBlot T48** or **BeeBlot**

for 20/30 teststrips:

**Dynablot Heat** or **Autoblot 3000H** or **B20**



## 6. Interpretation

Place processed and **dried** Teststrips into the designated fields of the Collector™ sheet.

Align them to the schematic drawing using the red marker line (top) and the green line (bottom) and fix them with adhesive tape.



# Optional: Automated recording, interpreting and electronic archiving of results

Scan Collector™ sheet and use StripAssay® Evaluator software for automated interpretation of band patterns.



**ViennaLab®** Online Calculator v2.17

INFO | ASSAYS | EVALUATE | MAIN WEBSITE

**β-Globin StripAssay® MED**  
Assay for the molecular analysis of frequent β-globin gene mutations  
REF-4-130, Revision 03/2010

[Modify Strip](#) [Change Assay](#)

SUMMARY

**IVS 1.1 [G>A] (c.92+1G>A) heterozygous**

DETAILS

— β-globin:

- -101 [C>T] (c.-151C>T) normal
- -87 [C>G] (c.-137C>G) normal
- -30 [T>A] (c.-80T>A) normal
- codon 5 [C>T] (c.17\_18delCT) normal
- codon 6 [G>A] H6C (c.19G>A) normal
- codon 6 [A>T] H6S (c.20A>T) normal
- codon 6 [A] (c.20delA) normal
- codon 8 [A>A] (c.25\_26delAA) normal
- codon 8/9 [T>G] (c.27\_28insG) normal
- codon 15 [TGG>TGA] (c.48G>A) normal
- codon 27 [G>T] Knossos (c.85G>T) normal
- IVS 1.1 [G>A] (c.92+1G>A) heterozygous
- IVS 1.5 [G>C] (c.92+5G>C) normal
- IVS 1.6 [T>C] (c.92+6T>C) normal
- IVS 1.110 [G>A] (c.93-21G>A) normal
- IVS 1.116 [T>G] (c.93-15T>G) normal
- IVS 1.130 [G>C] (c.93-1G>C) normal
- codon 39 [C>T] (c.118C>T) normal
- codon 44 [C] (c.135delC) normal
- IVS 2.1 [G>A] (c.315+1G>A) normal
- IVS 2.745 [C>G] (c.316-106C>G) normal
- IVS 2.848 [C>A] (c.316-3C>A) normal



**B|G| ViennaLab®**

ViennaLab Diagnostics GmbH  
Gaudenzdorfer Guertel 43-45  
A-1120 Vienna, Austria

phone: +43 1 8120156  
techhelp@viennalab.com  
www.viennalab.com

ISO 13485